

Equal Component Impression Material Comparison Chart

EQUAL COMPONENT SYSTEMS (2-PART)							
Material	Item Number	Color	Manufacturer	Number of Impressions	Viscosity Setting	Time	Shore-A Value
Gold Velvet II	GV-10	Yellow	JKR Labs	10	Low	6 min	40
Gold Velvet II	GV-40	Yellow	JKR Labs	40	Low	6 min	40
Gold Velvet II	GV-80	Yellow	JKR Labs	80	Low	6 min	40
Dur-A-Sil Equal	RB55	Blue	ipi	50-55	Medium	5-6 min	25
Vinyl Polysiloxane	VS-100/60	Green	Precision Labs	60	Medium	5-7 min	30
Silhouette	SILT	Blue	Siemens	55	Medium-High	6-8 min	26 +/- 2
Audiologist's Choice	01601	Teal	Oaktree Products	55	High	4m +/- 30s	39-40
Otoform A/K	OTO80	Green	Dreve	80	High	4m +/- 30s	35 +/- 2
Pink Silicast	245A	Pink	Westone	45	High	3-5 min	24-30
Pink Silicast	270A	Pink	Westone	70	High	3-5 min	24-30